

For Research Use Only

CENPT Polyclonal antibody, PBS Only

Catalog Number: 33475-1-PBS



Basic Information

Catalog Number:

33475-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38222

GenBank Accession Number:

BC015202

GeneID (NCBI):

80152

UNIPROT ID:

Q96BT3

Full Name:

centromere protein T

Calculated MW:

561 aa, 60 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33475-1-AP (CENPT antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33475-1-PBS in a different storage buffer formulation.



IP result of anti-CENPT (IP:33475-1-AP, 4ug; Detection:33475-1-AP 1:3000) with HEK-293T cells lysate 1840 ug. This data was developed using the same antibody clone with 33475-1-PBS in a different storage buffer formulation.