

For Research Use Only

Angiopoietin 4 Polyclonal antibody

Catalog Number:33510-1-AP



Basic Information

Catalog Number:

33510-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36152

GenBank Accession Number:

BC111976

GeneID (NCBI):

51378

UNIPROT ID:

Q9Y264

Full Name:

angiopoietin 4

Observed MW:

57 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HepG2 cells, L02 cells, NCI-H1299 cells, RAW 264.7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33510-1-AP (ANGPT4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.