

For Research Use Only

FAM108C1 Polyclonal antibody, PBS Only

Catalog Number:33561-1-PBS



Basic Information

Catalog Number:

33561-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37831

GenBank Accession Number:

BC059401

GeneID (NCBI):

58489

UNIPROT ID:

Q6PCB6

Full Name:

family with sequence similarity 108, member C1

Calculated MW:

329 aa, 36 kDa

Observed MW:

36 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33561-1-AP (FAM108C1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33561-1-PBS in a different storage buffer formulation.