

For Research Use Only

ZNF334 Polyclonal antibody

Catalog Number:33567-1-AP



Basic Information

Catalog Number:

33567-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36359

GenBank Accession Number:

BC024177

GeneID (NCBI):

55713

UNIPROT ID:

Q9HCZ1

Full Name:

zinc finger protein 334

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human PBMCs,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33567-1-AP (ZNF334 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.