

For Research Use Only

ZNF334 Polyclonal antibody, PBS Only

Catalog Number:33567-1-PBS



Basic Information

Catalog Number:

33567-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36359

GenBank Accession Number:

BC024177

GeneID (NCBI):

55713

UNIPROT ID:

Q9HCZ1

Full Name:

zinc finger protein 334

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33567-1-AP (ZNF334 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33567-1-PBS in a different storage buffer formulation.