

For Research Use Only

C9orf114 Polyclonal antibody, PBS Only

Catalog Number: 33571-1-PBS



Basic Information

Catalog Number: 33571-1-PBS	GenBank Accession Number: BC063644	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 51490	
Source: Rabbit	UNIPROT ID: Q5T280	
Isotype: IgG	Full Name: chromosome 9 open reading frame 114	
Immunogen Catalog Number: AG38014	Observed MW: 40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

SPOUT1 (SPOUT domain-containing methyltransferase 1), also known as CENP32 or C9orf114, is a methyltransferase containing a SPOUT domain. This protein is highly expressed in the nucleolus and is crucial for ribosome function. SPOUT1 has been identified as an essential protein for cell survival in both Drosophila and humans. Moreover, mutations in SPOUT1 are associated with developmental epileptic encephalopathy (DEE) in humans, characterized by infantile seizures and intellectual disability.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33571-1-AP (C9orf114 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33571-1-PBS in a different storage buffer formulation.