

For Research Use Only

CNTN5 Polyclonal antibody

Catalog Number:33610-1-AP



Basic Information

Catalog Number:

33610-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38223

GenBank Accession Number:

BC130619

GeneID (NCBI):

53942

UNIPROT ID:

O94779

Full Name:

contactin 5

Calculated MW:

121 kDa

Observed MW:

121-150 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33610-1-AP (CNTN5 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.