

For Research Use Only

Histone H2B Polyclonal antibody

Catalog Number:33612-1-AP



Basic Information

Catalog Number: 33612-1-AP	GenBank Accession Number: NM_021063.3	Purification Method: Antigen affinity Purification
Source: Rabbit	GeneID (NCBI): 3017	Recommended Dilutions: WB: 1:1000-1:4000 IHC: 1:200-1:800
Isotype: IgG	UNIPROT ID: P58876	
Immunogen Catalog Number: AG38954	Full Name: histone cluster 1, H2bd	
	Calculated MW: 14kDa, 126aa	
	Observed MW: 16 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Species Specificity:
human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HEK-293 cells, HSC-T6 cells, HeLa cells, NIH/3T3 cells

IHC : mouse brain tissue,

Background Information

Histone H2B is a core component of the nucleosome, the fundamental unit of chromatin. It forms a dimer with Histone H2A, and two H2A-H2B dimers combine with a (H3-H4)₂ tetramer to create the histone octamer. The primary function of this structure is to package DNA into a compact and manageable form. However, H2B is also subject to various post-translational modifications (PTMs), such as ubiquitination, which serve as epigenetic marks to regulate gene expression by altering chromatin accessibility.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33612-1-AP (Histone H2B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 33612-1-AP (Histone H2B antibody) at dilution of 1:400 (under 40x Lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).