

For Research Use Only

TRABD Polyclonal antibody, PBS Only

Catalog Number: 33633-1-PBS



Basic Information

Catalog Number:

33633-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40480

GenBank Accession Number:

BC093029

GeneID (NCBI):

80305

UNIPROT ID:

Q9H413

Full Name:

TraB domain containing

Observed MW:

37 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

TraB domain-containing protein (TRABD) is also named TTG2. TRABD is a Tiki/TraB family protein that is highly expressed in the brains of patients with AD (PMID: 33446646). TRABD as a novel outer mitochondrial membrane protein. Depletion of TRABD in cells severely impairs mitochondrial respiration and ATP production, inhibits cell growth, increases reactive oxygen species levels (PMID: 40778267).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33633-1-AP (TRABD antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33633-1-PBS in a different storage buffer formulation.



IP result of anti-TRABD (IP:33633-1-AP, 4ug; Detection:33633-1-AP 1:1000) with HeLa cells lysate 920 ug. This data was developed using the same antibody clone with 33633-1-PBS in a different storage buffer formulation.