

For Research Use Only

# C8orf55 Polyclonal antibody

Catalog Number:33642-1-AP



## Basic Information

Catalog Number:

33642-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38944

GenBank Accession Number:

NM\_016647.2

GeneID (NCBI):

51337

UNIPROT ID:

Q8WUY1

Full Name:

chromosome 8 open reading frame 55

Calculated MW:

24kDa, 208aa

Observed MW:

23 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MCF-10A cells, MCF-7 cells, MDA-MB-231 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33642-1-AP (C8orf55 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.