

For Research Use Only

C8orf55 Polyclonal antibody, PBS Only

Catalog Number:33642-1-PBS



Basic Information

Catalog Number:

33642-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38944

GenBank Accession Number:

NM_016647.2

GeneID (NCBI):

51337

UNIPROT ID:

Q8WUY1

Full Name:

chromosome 8 open reading frame 55

Calculated MW:

24kDa, 208aa

Observed MW:

23 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33642-1-AP (C8orf55 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33642-1-PBS in a different storage buffer formulation.