

For Research Use Only

MUC15 Polyclonal antibody

Catalog Number:33669-1-AP



Basic Information

Catalog Number: 33669-1-AP	GenBank Accession Number: NM_145650	Purification Method: Antigen affinity Purification
Source: Rabbit	GeneID (NCBI): 143662	Recommended Dilutions: WB: 1:500-1:1000
Isotype: IgG	UNIPROT ID: Q8N387-1	
Immunogen Catalog Number: EG5625	Full Name: mucin 15, cell surface associated	
	Calculated MW: 36 kDa	
	Observed MW: 85 kD, 110 kD	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue, PNGF treated human placenta tissue
Species Specificity: human	

Background Information

Mucin 15(MUC15) is a member of the transmembrane glycoprotein mucins that are abundantly expressed in the epithelial cells of mammals. MUC15 is involved in cell adhesion and signal transduction, and it plays a role in promoting or inhibiting tumor progression. (PMID: 12047385; 17720698;21615302)

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33669-1-AP (MUC15 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.