

For Research Use Only

TMEM260 Polyclonal antibody

Catalog Number:33682-1-AP



Basic Information

Catalog Number:

33682-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38950

GenBank Accession Number:

NM_017799.3

GeneID (NCBI):

54916

UNIPROT ID:

Q9NX78

Full Name:

chromosome 14 open reading frame
101

Calculated MW:

80 kDa, 707aa

Observed MW:

80 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, HuH-7 cells, L02 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33682-1-AP (TMEM260 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.