

For Research Use Only

# CNTN3 Polyclonal antibody

Catalog Number:33720-1-AP



## Basic Information

Catalog Number:

33720-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5928

GenBank Accession Number:

NM\_020872.2

GeneID (NCBI):

5067

UNIPROT ID:

Q9P232

Full Name:

contactin 3 (plasmacytoma associated)

Calculated MW:

113 kDa

Observed MW:

140 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse brain tissue, mouse cerebellum tissue, rat brain tissue, rat cerebellum tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33720-1-AP (CNTN3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.