

For Research Use Only

Contactin-6 Polyclonal antibody

Catalog Number:33731-1-AP



Basic Information

Catalog Number:

33731-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5753

GenBank Accession Number:

NM_001289080.2

GeneID (NCBI):

27255

UNIPROT ID:

Q9UQ52

Full Name:

contactin 6

Calculated MW:

114 kDa

Observed MW:

114 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse cerebellum tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse cerebellum tissue lysate was subjected to SDS PAGE followed by western blot with 33731-1-AP (Contactin-6 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.