

For Research Use Only

ZCCHC24 Polyclonal antibody

Catalog Number:33741-1-AP



Basic Information

Catalog Number:

33741-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG39605

GenBank Accession Number:

NM_153367.4

GeneID (NCBI):

219654

UNIPROT ID:

Q8N2G6

Full Name:

zinc finger, CCHC domain containing
24

Calculated MW:

27kDa, 241aa

Observed MW:

27 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : LNCaP cells, PC-3 cells, U2OS cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33741-1-AP (ZCCHC24 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.