

For Research Use Only

SLC35A5 Polyclonal antibody

Catalog Number: 33755-1-AP



Basic Information

Catalog Number: 33755-1-AP	GenBank Accession Number: BC005207	Purification Method: Antigen affinity Purification
Source: Rabbit	GeneID (NCBI): 55032	Recommended Dilutions: WB: 1:1000-1:4000
Isotype: IgG	UNIPROT ID: Q9BS91	
Immunogen Catalog Number: AG8748	Full Name: solute carrier family 35, member A5	
	Calculated MW: 424 aa, 49 kDa	
	Observed MW: 45-50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : DU 145 cells, PC-3 cells
Species Specificity: human	

Background Information

Solute carrier family 35 member A5 (SLC35A5) is a member of the SLC35A protein subfamily comprising nucleotide sugar transporters. SLC35A5 is the sole member of the SLC35 family that contains several diacidic motifs within its cytosolically exposed C-terminus, some of which might be capable of UDP-sugar binding (PMID: 30641943).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33755-1-AP (SLC35A5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.