

For Research Use Only

Siglec-H Polyclonal antibody

Catalog Number: 33801-1-AP



Basic Information

Catalog Number:

33801-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6216

GenBank Accession Number:

/

GeneID (NCBI):

233274

UNIPROT ID:

Q3Y597

Full Name:

sialic acid binding Ig-like lectin H

Calculated MW:

34 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : mouse kidney tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse kidney tissue lysates were subjected to SDS PAGE followed by western blot with 33801-1-AP (Siglec-H antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.