

For Research Use Only

C10orf35 Polyclonal antibody

Catalog Number:33819-1-AP



Basic Information

Catalog Number:

33819-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36807

GenBank Accession Number:

BC013587

GeneID (NCBI):

219738

UNIPROT ID:

Q96D05

Full Name:

chromosome 10 open reading frame

35

Observed MW:

15 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MG-63 cells, U-251 cells, U2OS cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33819-1-AP (C10orf35 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.