

For Research Use Only

PADI3 Polyclonal antibody

Catalog Number:33848-1-AP



Basic Information

Catalog Number:

33848-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35821

GenBank Accession Number:

NM_016233.2

GeneID (NCBI):

51702

UNIPROT ID:

Q9ULW8

Full Name:

peptidyl arginine deiminase, type III

Calculated MW:

75kDa, 664aa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse kidney tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33848-1-AP (PADI3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.