

For Research Use Only

ACBD7 Polyclonal antibody, PBS Only

Catalog Number:33855-1-PBS



Basic Information

Catalog Number:

33855-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40829

GenBank Accession Number:

BC029526

GeneID (NCBI):

414149

UNIPROT ID:

Q8N6N7

Full Name:

acyl-Coenzyme A binding domain containing 7

Calculated MW:

88 aa, 10 kDa

Observed MW:

10 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data