

For Research Use Only

# PID1 Polyclonal antibody, PBS Only

Catalog Number:33859-1-PBS



## Basic Information

Catalog Number:

33859-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40871

GenBank Accession Number:

BC040164

GeneID (NCBI):

55022

UNIPROT ID:

Q7Z2X4

Full Name:

phosphotyrosine interaction domain containing 1

Calculated MW:

168 aa, 19 kDa

Observed MW:

25 kDa, 28 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

## Background Information

PID1 is a cytoplasmic PTB-domain adaptor that negatively regulates insulin and IGF1 signaling, modulating AKT activity, glucose uptake, and mitochondrial function. Enriched in muscle, adipose tissue, and liver, PID1 promotes insulin resistance, impacts lipid handling, and contributes to metabolic disease. Its dysregulation also influences tumor metabolism and stress responses, underscoring PID1 as an important node linking growth factor signaling, energy homeostasis, and pathophysiology.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33859-1-AP (PID1 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33859-1-PBS in a different storage buffer formulation.