

For Research Use Only

ROPN1B Polyclonal antibody

Catalog Number:33925-1-AP



Basic Information

Catalog Number:

33925-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40678

GenBank Accession Number:

BC141849

GeneID (NCBI):

152015

UNIPROT ID:

Q9BZX4

Full Name:

ropporin, rhophilin associated protein
1B

Calculated MW:

24 kDa

Observed MW:

22 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human testis tissue, mouse testis tissue, rat testis
tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33925-1-AP (ROPN1B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.