

For Research Use Only

# CENPS-CORT Polyclonal antibody, PBS Only

Catalog Number:33963-1-PBS



## Basic Information

Catalog Number:

33963-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG40952

GenBank Accession Number:

GeneID (NCBI):

100526739

UNIPROT ID:

Q8N2Z9

Full Name:

CENPS-CORT readthrough

Calculated MW:

15kDa

Observed MW:

19 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

CENPS-CORT is a centromeric protein complex composed of core subunits such as CENPS and CENPT, playing a critical role in chromosome biology. This complex localizes to the centromere region and directly participates in the assembly and maintenance of the centromere-kinetochore structure, ensuring accurate chromosome segregation during both mitosis and meiosis. Its functions primarily include promoting nucleosome assembly at the centromere, maintaining the three-dimensional structural integrity of the centromere, and ensuring the fidelity of chromosome separation by regulating microtubule-kinetochore attachments. Studies have shown that dysfunction of the CENPS-CORT complex disrupts the accuracy of chromosome segregation, leading to chromosomal numerical abnormalities such as aneuploidy, which is closely associated with tumorigenesis and reproductive disorders. It is therefore an important subject in research on cell cycle regulation and genomic stability.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33963-1-AP (CENPS-CORT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33963-1-PBS in a different storage buffer formulation.