

For Research Use Only

Mouse MD-1 Polyclonal antibody, PBS Only

Catalog Number: 34213-1-PBS



Basic Information

Catalog Number:

34213-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6369

GenBank Accession Number:

NM_010745.2

GeneID (NCBI):

17084

UNIPROT ID:

O88188

Full Name:

Lymphocyte antigen 86

Calculated MW:

18 kDa

Observed MW:

24 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

mouse, rat

Background Information

Lymphocyte antigen 86 (LY86), also known as myeloid differentiation protein 1 (MD1), encodes a secreted glycoprotein that participates in inflammation, humoral immune response, apoptosis, and signal transduction. LY86 is widely expressed in cardiac tissues and predominantly by B cells and macrophages. (PMID: 38822281)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Rat spleen tissue were subjected to SDS PAGE followed by western blot with 34213-1-AP (MD-1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 34213-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded rat spleen tissue slide using 34213-1-AP (MD-1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 34213-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded rat spleen tissue slide using 34213-1-AP (MD-1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 34213-1-PBS in a different storage buffer formulation.