

For Research Use Only

GFP Monoclonal antibody (3H9)



www.ptgcn.com

Catalog Number: 3h9-2

Basic Information

Catalog Number:
3h9-2

Applications:
WB

Host/IsoType:
Rat/IgG2a

Conjugate:
Unconjugated

Type:
Primary Antibody

RRID:
AB_10773374

Description

Rat monoclonal antibody [3H9] to Green Fluorescent Proteins

Specificity/Target

The antibody recognizes green fluorescent proteins derived from *Aequorea victoria* GFP, such as eGFP, wtGFP, YFP, or CFP.

Recommended Dilution

Western Blot: 1:1,000

Storage

Storage:

Shipped on ice. Store at -20°C/-4°F. Stable for 1 year after shipment.

Storage Buffer:

PBS, 50% glycerol, preservative: 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Dilution series of EGFP in HEK 293T cell lysate
Primary antibody: [3H9] 1:1,000 Secondary
antibody: goat anti-rat Alexa Fluor® 488 1:1,000