

For Research Use Only

ABCE1 Polyclonal antibody

Catalog Number: 51008-2-AP



Basic Information

Catalog Number:

51008-2-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0454

GenBank Accession Number:

BC016283

GeneID (NCBI):

6059

UNIPROT ID:

P61221

Full Name:

ATP-binding cassette, sub-family E (OABP), member 1

Calculated MW:

599 aa, 67 kDa

Observed MW:

67 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

