

For Research Use Only

# Evi2a Polyclonal antibody

Catalog Number: 51014-2-AP



## Basic Information

Catalog Number:

51014-2-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0478

GenBank Accession Number:

BC038124

GeneID (NCBI):

14017

UNIPROT ID:

P20934

Full Name:

ecotropic viral integration site 2a

Calculated MW:

24 kDa

Observed MW:

26 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : mouse liver tissue,

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 51014-2-AP (Evi2a antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.