

For Research Use Only

# Harmonin Polyclonal antibody

Catalog Number: 51028-2-AP

1 Publications



## Basic Information

Catalog Number:

51028-2-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0513

GenBank Accession Number:

BC016057

GeneID (NCBI):

10083

UNIPROT ID:

Q9Y6N9

Full Name:

Usher syndrome 1C (autosomal recessive, severe)

Calculated MW:

552 aa, 62 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

IHC

Species Specificity:

human, mouse

Cited Species:

mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Grimsley-Myers Cynthia M CM	20016102	J Neurosci	IHC

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

