

For Research Use Only

# SSTR5 Polyclonal antibody

Catalog Number: 51029-2-AP



## Basic Information

Catalog Number:

51029-2-AP

Size:

347 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0519

GenBank Accession Number:

BC146576

GeneID (NCBI):

6755

UNIPROT ID:

P35346

Full Name:

somatostatin receptor 5

Calculated MW:

364 aa, 39 kDa

Observed MW:

50 kDa, 65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human placenta tissue,

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 51029-2-AP (SSTR5 antibody) at dilution of 1:0 incubated at room temperature for 1.5 hours.