

For Research Use Only

PROX1 Polyclonal antibody

Catalog Number: 51043-1-AP

Featured Product

10 Publications



Basic Information

Catalog Number:

51043-1-AP

Size:

150 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0457

GenBank Accession Number:

BC024201

GeneID (NCBI):

5629

UNIPROT ID:

Q92786

Full Name:

prospero homeobox 1

Calculated MW:

737 aa, 83 kDa

Observed MW:

90 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

ChIP, CoIP, IF, IHC, WB

Species Specificity:

human, mouse, rat

Cited Species:

hamster, human, mouse

Background Information

PROX1, Prospero homeobox protein 1, belongs to the Prospero homeobox family and contain a Prospero-type homeobox DNA-binding domain. PROX1 plays a fundamental role in early development of CNS and regulated the expression of gene and development of postmitotic undifferentiated young neuron. Also PROX1 is a specific and required regulator of the development of the lymphatic system. This is a rabbit polyclonal antibody raised against part of C-terminal chain of human PROX1.

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-----------------|-----------|------------------|-------------|
| Sean M Armour | 28916805 | Nat Commun | CoIP, chIP |
| Xu Shang-Zhi SZ | 23675176 | Int J Biomed Sci | WB |
| Yuki Kawai | 24680387 | J Surg Res | IF |

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data