

For Research Use Only

PNRC1 Polyclonal antibody

Catalog Number: 51052-1-AP

1 Publications



Basic Information

Catalog Number:

51052-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0100

GenBank Accession Number:

BC018112

GeneID (NCBI):

10957

UNIPROT ID:

Q12796

Full Name:

proline-rich nuclear receptor
coactivator 1

Calculated MW:

35 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB: Raji cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xia Shen	39046443	Elife	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Raji cells were subjected to SDS PAGE followed by western blot with 51052-1-AP (PNRC1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.