

For Research Use Only

Hepsin Polyclonal antibody

Catalog Number: 51060-1-Ig



Basic Information

Catalog Number:

51060-1-Ig

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC138809

GeneID (NCBI):

15451

UNIPROT ID:

O35453

Full Name:

hepsin

Calculated MW:

47 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human liver tissue, mouse liver tissue

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 51060-1-Ig (Hepsin antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.