

For Research Use Only

# Hepsin Polyclonal antibody

Catalog Number: 51060-1-Ig



## Basic Information

Catalog Number:

51060-1-Ig

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC138809

GeneID (NCBI):

15451

UNIPROT ID:

O35453

Full Name:

hepsin

Calculated MW:

47 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human liver tissue, mouse liver tissue

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human liver tissue were subjected to SDS PAGE followed by western blot with 51060-1-Ig (Hepsin antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.