

For Research Use Only

# JAK1 Polyclonal antibody

Catalog Number: 51079-1-AP

6 Publications



## Basic Information

Catalog Number:

51079-1-AP

Concentration:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC132729

GeneID (NCBI):

3716

UNIPROT ID:

P23458

Full Name:

Janus kinase 1 (a protein tyrosine kinase)

Calculated MW:

1154 aa, 133 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, rat

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jiangang Zheng	36572890	BMC Vet Res	WB
Long Zhou	39669585	Front Immunol	WB
Muzi Li	39209829	Cell Death Discov	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

