

For Research Use Only

MAP4K4 Polyclonal antibody

Catalog Number: 51088-1-AP **1 Publications**



Basic Information

Catalog Number:

51088-1-AP

Size:

290 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC156887

GeneID (NCBI):

9448

UNIPROT ID:

O95819

Full Name:

mitogen-activated protein kinase
kinase kinase kinase 4

Calculated MW:

142 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Lifeng Ding	37647682	Int Immunopharmacol	WB,IHC

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data