

For Research Use Only

PLCB3 Polyclonal antibody

Catalog Number: 51139-1-AP



Basic Information

Catalog Number:

51139-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC142681

GeneID (NCBI):

5331

UNIPROT ID:

Q01970

Full Name:

phospholipase C, beta 3
(phosphatidylinositol-specific)

Calculated MW:

1234 aa, 139 kDa

Observed MW:

139 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HEK-293T cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 51139-1-AP (PLCB3 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.