

For Research Use Only

# SOX21 Polyclonal antibody

Catalog Number: 55052-1-AP

3 Publications



## Basic Information

Catalog Number:

55052-1-AP

Size:

1300 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_007084

GeneID (NCBI):

11166

UNIPROT ID:

Q9Y651

Full Name:

SRX (sex determining region Y)-box  
21

Calculated MW:

29 kDa

Observed MW:

25-30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, goat

Positive Controls:

WB : mouse embryo tissue, Jurkat cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Hailey E Brighton	29180473	Cancer Res	WB
Sen Ma	32806500	Animals (Basel)	WB
Qiaochang Cao	38992573	Life Sci	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse embryo tissue were subjected to SDS PAGE followed by western blot with 55052-1-AP (SOX21 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.