

For Research Use Only

# MTCO1 Polyclonal antibody

Catalog Number: 55071-1-AP

4 Publications



## Basic Information

Catalog Number:

55071-1-AP

Size:

290 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

YP\_003024028

GeneID (NCBI):

4512

Full Name:

cytochrome c oxidase I

Calculated MW:

57 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Jing Zhang	26492917	EMBO J	WB
Hao-Peng Lin	27165101	Biomed Microdevices	WB
Xiaoting Liu	35876654	EMBO Rep	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

