

For Research Use Only

# CSPG5 Polyclonal antibody

Catalog Number: 55086-1-AP



## Basic Information

Catalog Number:

55086-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_006574

GeneID (NCBI):

10675

UNIPROT ID:

O95196

Full Name:

chondroitin sulfate proteoglycan 5  
(neuroglycan C)

Calculated MW:

60 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

CSPG5, also named as CALEB and NGC, may function as a growth and differentiation factor involved in neuritogenesis. It may induce ERBB3 activation. The antibody is specific to CSPG5.

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

