

For Research Use Only

GRIN2D Polyclonal antibody

Catalog Number: 55094-1-AP



Basic Information

Catalog Number:

55094-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000836

GeneID (NCBI):

2906

UNIPROT ID:

O15399

Full Name:

glutamate receptor, ionotropic, N-methyl D-aspartate 2D

Calculated MW:

144 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

