

For Research Use Only

# CCNB3 Polyclonal antibody

Catalog Number: 55095-1-AP



## Basic Information

Catalog Number:

55095-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_033031

GeneID (NCBI):

85417

UNIPROT ID:

Q8WWL7

Full Name:

cyclin B3

Calculated MW:

158 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Background Information

CCNB3, also named as CYCB3, belongs to the cyclin family and Cyclin AB subfamily. Cyclins are positive regulatory subunits of the cyclin-dependent kinases (CDKs), and thereby play an essential role in the control of the cell cycle, notably via their destruction during cell division. Its tissue specificity suggest that it may be required during early meiotic prophase I. The antibody is specific to CCNB3.

## Storage

Storage:

Store at -20°C.

Storage Buffer:

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data