

For Research Use Only

# IL-11 Polyclonal antibody

Catalog Number: 55169-1-AP

Featured Product

9 Publications



## Basic Information

Catalog Number:

55169-1-AP

Size:

650 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000641

GeneID (NCBI):

3589

UNIPROT ID:

P20809

Full Name:

interleukin 11

Calculated MW:

21 kDa

Observed MW:

25 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, monkey

Positive Controls:

WB : human plasma tissue,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Xin Bai	35731866	Sci Adv	WB, IHC
Yuchen Li	31864417	Vet Res	WB
Fengyan Huang	31841894	Ecotoxicol Environ Saf	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with 55169-1-AP (IL11 antibody) at dilution of 1:1000 incubated at 4 degree celsius over night.