

For Research Use Only

BNC2 Polyclonal antibody

Catalog Number: 55220-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

55220-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_017637

GeneID (NCBI):

54796

UNIPROT ID:

Q6ZN30

Full Name:

basonuclin 2

Calculated MW:

122 kDa

Observed MW:

122-140 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : HeLa cells, mouse uterus tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Marie Bobowski-Gerard	36088459	Nat Commun	WB,IP
Ninon Very	38830870	Cell Death Dis	IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 55220-1-AP (BNC2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.