

For Research Use Only

# BNC1 Polyclonal antibody

Catalog Number: 55233-1-AP



## Basic Information

Catalog Number:

55233-1-AP

Size:

413 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001717

GeneID (NCBI):

646

UNIPROT ID:

Q01954

Full Name:

basonuclin 1

Calculated MW:

111 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : HepG2 cells,

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using 55233-1-AP (BNC1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).