

For Research Use Only

MYD88 Polyclonal antibody

Catalog Number: 55250-1-AP **9 Publications**



Basic Information

Catalog Number:

55250-1-AP

Concentration:

140 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002468

GeneID (NCBI):

4615

UNIPROT ID:

Q99836

Full Name:

myeloid differentiation primary response gene (88)

Calculated MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, rat

Positive Controls:

WB : K-562 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jiali Shi	39892170	Int Immunopharmacol	WB
Qin Yang	39617843	Sci Rep	WB
Wenyan Yu	39520083	Adv Sci (Weinh)	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

