

For Research Use Only

BCL2L11, BIM Polyclonal antibody

Catalog Number: 55269-1-AP

1 Publications



Basic Information

Catalog Number:

55269-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_006538

GeneID (NCBI):

10018

UNIPROT ID:

O43521

Full Name:

BCL2-like 11 (apoptosis facilitator)

Calculated MW:

23 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Binhui Pan	38831322	BMC Med Genomics	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

