

For Research Use Only

DOT1L Polyclonal antibody

Catalog Number: 55315-1-AP



Basic Information

Catalog Number:

55315-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_032482

GeneID (NCBI):

84444

UNIPROT ID:

Q8TEK3

Full Name:

DOT1-like, histone H3
methyltransferase (S. cerevisiae)

Calculated MW:

185 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

