

For Research Use Only

MTIF2 Polyclonal antibody, PBS Only

Catalog Number:55325-1-PBS



Basic Information

Catalog Number:

55325-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002453

GeneID (NCBI):

4528

UNIPROT ID:

P46199

Full Name:

mitochondrial translational initiation factor 2

Calculated MW:

81 kDa

Observed MW:

81 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

MTIF2 belongs to the IF-2 family. It is one of the essential components for the initiation of protein synthesis. MTIF2 protects formylmethionyl-tRNA from spontaneous hydrolysis and promotes its binding to the 30S ribosomal subunits. It is also involved in the hydrolysis of GTP during the formation of the 70S ribosomal complex. This antibody is specific to MTIF2.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human heart tissue were subjected to SDS PAGE followed by western blot with 55325-1-AP (MTIF2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 55325-1-PBS in a different storage buffer formulation.