

For Research Use Only

QSER1 Polyclonal antibody, PBS Only

Catalog Number:55326-1-PBS



Basic Information

Catalog Number:

55326-1-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001076786

GeneID (NCBI):

79832

UNIPROT ID:

Q2KHR3

Full Name:

glutamine and serine rich 1

Calculated MW:

190 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IHC, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

QSER1, a recently characterized protein, whose activity has been tightly associated with the protection of DNA methylation sites in hESCs. QSER1 is widely up-regulated in multiple cancers, and its up-regulation correlates with poor clinic outcomes (PMID: 36371602). QSER1 depletion leads to NODAL de-repression and associates with increased RNA Polymerase II occupancy at YAP1bound enhancers, pointing to a QSER1 contribution maintaining low transcriptional activity on certain genomic regions (PMID: 33833093, 35210841).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 55326-1-AP (QSER1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 55326-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 55326-1-AP (QSER1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 55326-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using QSER1 antibody (55326-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 55326-1-PBS in a different storage buffer formulation.